

## PORK SAFETY

**Title:** Rapid detection of *Yersinia enterocolitica* in pigs using the TaqMan System - **NPB #98/160**

**Investigator:** Irene V. Wesley

**Institution:** USDA-ARS-National Animal Disease Center  
2300 Dayton Road  
Ames, Iowa 50010

**Date Received:** 10/15/1999

**ABSTRACT:** The TaqMan assay, or 5' nuclease PCR assay, is a second-generation PCR detection system, which is reportedly more sensitive than conventional PCR tests. We have developed and evaluated a 5' nuclease PCR assay for the detection of *Y. enterocolitica*. The assay targets the chromosomally encoded *ail* (adhesion invasion locus) gene.

Three primer/probe sets, (TM1, TM2, and TM3) amplifying different, yet overlapping, regions of *ail* were examined for specificity and sensitivity. The TM1 set displayed the highest specificity, accurately detecting each of the 26 *Y. enterocolitica* strains and none of the 21 non-*enterocolitica* strains. TM1 set detected ~0.5 pg ( $10^{-12}$  grams) of purified *Y. enterocolitica* DNA. The TM2 set was the most sensitive and detected ~ 0.25 pg of purified DNA. However, it failed to recognize 10 of the *Y. enterocolitica* strains used in this study. For TM3, sensitivities comparable to TM1 were achieved; cross-reaction with non-*enterocolitica* strains was not observed. However, TM3 did not identify all of the *Y. enterocolitica* strains tested.

The optimized TaqMan assay was compared with bacteriological culture methods and the first-generation multiplex PCR for the rapid detection of pathogenic *Y. enterocolitica* in market weight hogs (n=240) and pork products (n=650). *Y. enterocolitica* was not detected by bacteriological culture in any of the hog tissues tested (nine samples per hog) but was detected by multiplex PCR (2.0%) and TaqMan (45.6%) assays. In addition, ground pork (n=300 samples) and chitterlings (n=350) were screened for *Y. enterocolitica*. By standard culture, *Y. enterocolitica* was detected in chitterlings (8%), but not in ground pork (0%). By multiplex PCR, *Y. enterocolitica* was identified in ground pork (12%) and chitterlings (27%). In contrast, the highly specific TaqMan assay identified *Y. enterocolitica* in ground pork (52%) and chitterlings (79%).

The results of this study indicate that the TaqMan probes and primers for the *ail* gene (TM1) are more specific for pathogenic strains of *Y. enterocolitica* than either bacteriological culture or first-generation multiplex PCR assays.

*These research results were submitted in fulfillment of checkoff funded research projects. This report is published directly as submitted by the project's principal investigator. This report has not been peer reviewed*

**For more information contact:**

**National Pork Board, P.O. Box 9114, Des Moines, Iowa USA**

800-456-7675, Fax: 515-223-2646, E-Mail: [porkboard@porkboard.org](mailto:porkboard@porkboard.org), Web: <http://www.porkboard.org/>